


The First Open Oxygen System

OxyChin
An OxyArm™ Product



Delivers Freedom

OxyChin™

www.oxyarm.com

Patented Technology

The Science - Predictable Oxygen Delivery

The Innovation - Directional Diffuser Technology

The Results - Open Oxygen Systems

SAFE

- When emesis occurs
- Provides quick access for oral suctioning
- Ideal for conscious sedation & post op oxygenation

EFFICIENT

- OxyChin delivers 1 - 15 LPM
- Rate Trend Monitoring with O₂/CO₂ model

COMFORTABLE

- Reduces claustrophobia
- Open design allows communication between patient and caregiver
- Minimal contact oxygen delivery system

ORDER CODES:

OC-1025 - box of 25

OC-1025-8 - with tubing - box of 25

OC-2025 - Rate Trend Monitoring - box of 25

OC-2025-8 - with tubing - Rate Trend Monitoring - box of 25

Oxygen Flow Settings	Approximate O ₂ Concentrations
1 LPM	21% - 27%
2 LPM	28% - 31%
4 LPM	32% - 35%
6 LPM	36% - 39%
8 LPM	40% - 43%
10 LPM	42% - 47%
12-15 LPM	48% - 50%

Manufactured by

Southmedic

50 Alliance Blvd., Barrie Ontario, Canada L4M 5K3

705-726-9383

705-728-9537 FAX

www.southmedic.com

Email: smedic@southmedic.on.ca

ISO 13485

EDI 406000117

P0017A-Nov/05(Mar.06)